

### Sample requirements

Library type	Material	Library prep kit	Quantity	Volume approx.	Quality evaluation - user
Whole-genome seq	Genomic DNA	KAPA Hyperprep	250 ng/ul*	20 µl	DNA needs to be run on a 1% agarose gel to confirm integrity.
ChIP-Seq	DNA ChIP-enriched	KAPA Hyperprep	all	20 µl	qPCR analysis is recommended to confirm enrichment.
Cut&Run	DNA Cut&RUN	KAPA Hyperprep	all		
librairie sh/CRISPR	PCR amplicon	Nextera	5-50 ng/µl	10 µl	Migrate on gel.
PCR-Seq	PCR amplicon	Nextera	5-50 ng/µl	10 µl	Migrate on gel.
mRNA-Seq	Total RNA	KAPA Hyperprep RNA	250 ng/ul*	20 µl	Bioanalyzer will be run at the platform.
rRNA-deplété RNA-Seq	Total RNA	FAST select + KAPA hyperprep RNA	250 ng/ul*	15 µl	Bioanalyzer will be run at the platform.
Small RNA seq	Total RNA	Qiaseq small RNA	250 ng/ul*	15 µl	Use purification columns that retain small fragments. Ex: MirVana

\* This is the ideal concentration. We can make libraries with less material, contact us for details.

Resuspend samples in H<sub>2</sub>O. Samples can be in 1.5ml tubes clearly identified.

### Samples Quality Requirements

For RNA, RIN>8 is strongly suggested to assess proper RNA quality (for plants and FFPE tissues, RIN can be lower). Libraries for slightly degraded RNA samples should be prepared using ribo-depletion protocol to avoid 3' biases with a polyA enrichment. Nanodrop 260/280 ratio should be around 1,8-2,1 and 260/230 ratio should be higher than 1,5. Values outside of these ranges indicate contamination.

For DNA, gel must reflect a high DNA quality. Nanodrop 260/280 ratio should be close to 1,8-2,0 and the 260/230 ratio should be higher than 2,0.

Note that the Nanodrop frequently over-estimates RNA/DNA concentration.

We will QC your samples and libraries using a combination of Nanodrop, Qubit and Bioanalyzer. Libraries are quantified by Qubit, run on a BioAnalyzer, diluted to 10nM and normalized with qPCR quantification.

Lanes 1, 3, 4 show high quality gDNA. Lane 2 shows degraded gDNA.

